# THE FLEX BELT®

INSTRUCTION MANUAL

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ABDOMINAL MUSCLE TONER

BY **bmr**  
bio-medical research

[www.theflexbelt.com](http://www.theflexbelt.com)

**CONTRA-INDICATIONS, WARNINGS & PRECAUTIONS**

**Section 1** lists the important information which must be read before using your toning belt.

Your **THE FLEX BELT** by BMR product is indicated for the improvement of abdominal muscle tone, for strengthening of abdominal muscles, and for the development of a firmer abdomen. Your **THE FLEX BELT** by BMR product is suitable for use by all healthy adults; however, as with other forms of exercise, some care is needed when using it. Always follow the guidelines below and read the instruction manual before use. Some of the points below are gender-specific.

**CONTRA-INDICATIONS:****Do not use if:**

- You have an electronic implant (e.g. cardiac pacemaker or defibrillator) - as your **THE FLEX BELT** by BMR product may interfere with the proper functioning of the implanted stimulator) or suffer from any other heart problem.
- You have any other form of abdominal implant. The safety has not been established for the use of electronic muscle stimulators over abdominal implants.
- You are pregnant. The safety of electronic muscle stimulation (EMS) over the pregnant uterus has not been established.
- You suffer from cancer. The effects of electronic stimulation on cancerous tissue are unknown.
- You are under medical supervision for cognitive dysfunction as you may not be able to comply with safety instructions.
- The unit is in close proximity (e.g., 1 yard) to shortwave or microwave therapy equipment or you are connected to high-frequency surgical equipment (e.g. diathermy), because of the risk of device malfunction due to electronic interference.
- Wearing the device necessitates placement over areas at which drugs/medicines are administered (short-term or long-term) by injection (e.g., hormone treatments).

**WARNINGS:**

- The long-term effects of chronic electronic stimulation are unknown.
- Unit should not be used while sleeping, because you cannot react to safety information.

**Wait before using your THE FLEX BELT product until:**

- At least 6 weeks after the birth of your baby. You must consult your doctor before use.
- One month after an IUD contraceptive device (e.g., coil) has been fitted. You must consult your doctor before use.
- At least 3 months after having a Caesarean section. You must consult your doctor before use.
- The heavy days of your period have finished because vigorous abdominal exercise is not recommended at this time.

**When applying the GelPads™ and belt, always remember to:**

- Place the GelPads and belt **ONLY** on the abdomen as indicated in this manual.
- Not place the GelPads or belt on the front or sides of the neck, across or through the heart (i.e., one GelPad on the front of the chest and one on the back), in the genital region nor on the head, because of the risk of stimulating inappropriate muscles and organs.
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Avoid any recent scars, broken or inflamed skin, areas of infection or susceptibility to acne, thrombosis or other vascular problems (e.g., varicose veins), or any parts of the body where feeling is limited.
- Avoid areas of injury or restricted movement (e.g., fractures or sprains).
- Avoid placing the pads over metal implants.

**PRECAUTIONS:**

- Caution should be used for patients with suspected or diagnosed epilepsy.
- Caution should be used in the following:
  - When there is a tendency to bleed internally following an injury or fracture.
  - Following recent surgical procedures when electronic stimulation may disrupt the healing process
  - Over areas of skin which lack normal sensation
- Your **THE FLEX BELT** by BMR product is not intended for medical use, for the treatment of any medical condition, nor for any permanent physical changes.

**Get your doctor's permission before using your THE FLEX BELT product if:**

- You wish to use your THE FLEX BELT by BMR product and have a bad back. You should also ensure the intensity is kept low.
- You have any serious illness or injury not mentioned in this guide.
- You have recently undergone a surgical procedure.
- You take insulin for diabetes.
- You want to use it on a young child.
- You suffer from muscle or joint problems.
- You use the unit as part of a rehabilitation program.

**Adverse reactions:**

- Skin irritations and burns have been reported with the use of powered muscle stimulators.
- On very rare occasions, first-time users of EMS have reported feeling light-headed or faint. We recommend that you use the product while seated until you become accustomed to the sensation.

**To reposition GelPads during a session:**

- Always pause the program currently running, unfasten the belt and then refasten it behind your back after the GelPads' positions have been adjusted.

**After strenuous exercise or exertion:**

- Always use a lower intensity to avoid muscle fatigue.

**Contact THE FLEX BELT or an authorized dealer if:**

- Your unit is not working correctly. Do not use in the meantime.
- You experience any irritation, skin reaction, hypersensitivity, or other adverse reaction. You should, however, note that some reddening of the skin under the belt is quite normal during and for a short time after an exercise session.

**Note:**

- An effective session should not cause undue discomfort.

**Important:**

- For first time users, muscle stimulation can be an unusual sensation. We recommend that you begin in a seated position with low stimulation intensity settings to familiarize yourself with the sensation before progressing to higher intensity settings.
- Keep your unit out of the reach of children.
- The studs and GelPads must not be connected to other objects.
- Do not use your unit at the same time as any other device which transfers an electrical current into the body (e.g., another muscle stimulator).
- Do not touch the GelPads or metal studs while the unit is switched on.
- Do not use while driving, operating machinery or cycling.
- The belt is for single person use only.
- Cease using your unit if you are feeling light headed or faint. Consult your doctor if this happens.
- Do not over exert yourself when using muscle stimulation. Any workout should be at a level comfortable for you.
- Do not use your unit if you are wearing a belly-button ring. Remove the ring before you begin a session.
- THE FLEX BELT will not accept responsibility if the guidelines and instructions supplied with this unit are not followed.
- Although compliant with applicable EMC requirements, this device may still interfere with more sensitive equipment, please move away or switch off.

**Note:** If you are in any doubt about using THE FLEX BELT by BMR product for any reason, please consult your doctor before use.

## WELCOME

Congratulations on purchasing your THE FLEX BELT by BMR. Please read your user manual fully before using this product.

THE FLEX BELT by BMR is the future in body toning and wellness, giving you the body shape you desire. This revolutionary new concept allows a single hand-held controller to be attached to a range of garments which tone different areas of the body, saving you money on your next THE FLEX BELT purchase.

THE FLEX BELT by BMR uses clinically proven EMS technology to exercise all of the muscles of your abdomen, not just those under the belt. A signal is sent to the abdominal muscles and as the signal strength increases, the muscles contract. You can use your THE FLEX BELT by BMR almost anywhere - while relaxing at home, or even working at the office.

Use your THE FLEX BELT by BMR 5 times a week and in just 8 weeks, you can expect to achieve stronger, firmer, more toned abdominal muscles.

For best results we recommend that you use your THE FLEX BELT by BMR in conjunction with a normal, healthy diet and exercise.

### USER MANUAL CONTENTS

Your user manual is divided into sections to make it easy to use.

- Section 1**      Contra-Indications, Warnings & Precautions
- Section 2**      Introducing the THE FLEX BELT by BMR
- Section 3**      Setting up your THE FLEX BELT by BMR
- Section 4**      Your first THE FLEX BELT by BMR session
- Section 5**      Getting results with your THE FLEX BELT by BMR
- Section 6**      Help guide
- Section 7**      How does muscle stimulation work?
- Section 8**      Technical specifications
- Section 9**      Product warranty

If you have any difficulty setting up your unit or have any other questions or concerns about THE FLEX BELT by BMR call THE FLEX BELT Careline at:

Phone US: (855) FLEXGEAR (855) 353 - 9432

International: (310) 362 - 0581

Email: [info@flexgear.com](mailto:info@flexgear.com)

[www.theflexbelt.com](http://www.theflexbelt.com)

## INTRODUCING THE FLEX BELT BY BMR

Section 2 lists the pack contents and explains each item before you set up and begin.

You should find the following parts inside your pack.

### 1. Belt:

Your belt has been ergonomically designed and is fastened around your waist, under your clothes, during use. Designed for both men & women, fits waist sizes 24" – 47" The belt has a pocket for holding your unit during use.

### 2. Unit:

The unit generates signals which are sent, via the belt and pads, to your muscles to make them contract. Your unit is rechargeable and takes approximately three hours to charge fully.

### 2a. Unit Connector:

This connects the unit to the toning belt. Make sure the unit is fully connected to the belt.

**IMPORTANT:** Do not disconnect the unit from the toning belt until it has been switched off.

### 3. Pack of three adhesive pads.

These adhesive pads are placed onto your belt before use. When you wrap the belt around your waist, these pads must be in contact with your skin so that signals can be sent to your muscles.

### 4. THE FLEX BELT Battery Charger:

Plug the charger into an electrical outlet and then connect your unit to the charger to charge your unit's battery.

### 5. Instruction Manual:

This is your detailed guide to using your THE FLEX BELT by BMR.

### 6. Quick Start Guide:

This is an additional guide that you can use to help get you started.

## UNIT CONTROLS

### 1. On/Off Button (⏻)

Press and hold this button for 2 seconds to turn your unit on or off. You may also pause the exercise by briefly pressing this button.

### 2. Increase Training Intensity (▲)

Press and hold these buttons to increase the intensity of the exercise. The right-hand button controls the right-hand side of your belt. The left-hand button controls the left-hand side.

### 3. Decrease Training Intensity (▼)

Press and hold these buttons to decrease the intensity of the exercise on either side of the belt.

### 4. Program Button (P)

Press this button to select the toning program you wish to use. There are 10 programs in total, which are explained later in this manual (page 8).

### 5. Information Button (i)

Press this button to see your training information, such as the maximum intensity for your current or last session and the number of sessions you have completed.

## UNIT DISPLAY

- (50) Training intensity from the left-hand side of the belt.
- (45) Training intensity from the right-hand side of the belt.
- 10:08 Counts down the time left in the current session in minutes and seconds. Also used to display error messages (see page 10).
- || Appears when the exercise has been paused (you must press the on/off button briefly to restart the exercise).
- ✂ Indicates that the sound is off.
- 🔑 Keylock is active - you cannot change the intensity levels or the program accidentally.
- ⚠ This symbol indicates a poor contact between the unit and the belt or between the pads and the skin (see page 12 for more details).
- 🔋 Shows the battery power remaining.
- ☐ Number of sessions completed.
- # Shows the program currently running (1 - 10).
- 4 This symbol appears during each contraction phase.
- This symbol appears during each contraction phase.
- (92 87) Maximum intensity for the current/ last session.

## SETTING UP THE FLEX BELT BY BMR

Section 3 explains how to set up THE FLEX BELT by BMR before using it for the first time.

**Step 1 - Charge the battery.** Insert your charger into an electrical outlet and then connect your unit to it (Fig. a). It may take approximately 2-3 hours to fully charge the battery. Your unit is fully charged when all three sections of the battery icon are full. The battery icon will flash when the battery is low and needs to be recharged.

**NOTE:** Do not leave your unit connected to the charger when the battery is fully charged.

Fig. a

**Step 2 - Connect the unit to the belt.** The unit and belt must be connected for THE FLEX BELT by BMR to work. Place the belt on a flat surface with the gold stitching facing you. 'Click' the unit connector into the belt connector as shown in Fig. b.


**IMPORTANT:** Make sure the unit and belt are connected correctly or your THE FLEX BELT by BMR will not work - this is a safety feature.

Fig. b

**Step 3 - Remove the covers from the patterned side of the pads and place the pads on the belt.** Turn the belt over so that the 3 silver studs are facing upwards. Remove the three adhesive pads from their pack. There is 1 large square pad and 2 smaller oval pads. One side of each pad has a grid pattern on it while the other side is plain black.


(i) Remove the covers from the **patterned side** of the large pad.

(ii) Place the **patterned side** of this pad over the middle stud and between the white lines.

(iii) Now remove the covers from the **patterned side** of the two smaller pads.

(iv) Place the **patterned side** of the smaller pads firmly onto the belt, either side of the large pad.


(v) The location of the smaller pads depends on your waist size. In the example below, the smaller pads are placed nearest to the large, central pad. This is for a user whose waist size is 24"-31". Use the sizing guidelines below to help you choose your correct pad location.

### BELT SIZES

L - For waist sizes 41"-47" / 104-119cm

M - For waist sizes 32"-40" / 81-102cm


S - For waist sizes 24"-31" / 61-79cm

### IMPORTANT:

Before you progress, make sure that all three silver studs are completely covered by the pads. Press the edges of all three pads firmly onto the belt before use.



**Step 4 - Remove the covers from the black side of all three pads. Do not throw these covers away as you will need to put them back onto the pads at the end of your session.**

**NOTE:** You should only use THE FLEX BELT by BMR pads with your THE FLEX BELT by BMR unit.

**Step 5 - Position the belt on your body.**

Wrap the belt around your waist so that the pads are pressed firmly against your skin, positioning the large square pad over your navel (belly button). Stretch the ends of the belt around your waist (Fig. a) until the small pads are between your hip bone and ribs on either side of your body (Fig. b). Fasten the belt tightly around your waist, but not so tight as to cause discomfort. The pads are water-based, so you may find they are cool when placed on your skin.

You are now ready to begin using THE FLEX BELT by BMR.

## YOUR FIRST TRAINING SESSION

Section 4 explains how to use THE FLEX BELT by BMR.

Muscle toning can be an unusual sensation, but a pleasant one. It may tickle to begin with. We recommend that you use the product while seated until you become accustomed to the sensation. This progresses to a smooth muscle contraction as the training intensity increases. To get started, follow these simple steps:

**Step 1:** Press and hold the on/off button (⏻) for two seconds to switch your unit on (Fig. a).

**Step 2:** Select the program you wish to use by pressing the program button (P - Fig. b). Then to start the program, press and hold the increase intensity buttons (▲ - Fig. c) until you feel your muscles contracting. Always choose an intensity level at which you feel a strong but comfortable contraction. This will vary from person to person. You should try to reach an intensity level of 15 or higher in your first session. The increasing intensity level is displayed on the screen.

**Step 3:** Continue increasing the intensity throughout the session if possible. You will feel the muscle contractions getting stronger as the intensity increases. Remember the harder you work your muscles the better the results you will see, but always reduce the intensity if you feel any discomfort. Your unit's intensity range is 0-150.

Fig. a

Fig. b

Fig. c

### PAUSING A SESSION

If you wish to pause a session before it is finished, simply press the on/off button briefly. The display will show that the unit is paused (||). To resume the session, briefly press the on/off button again. The display will return to its normal mode.

### INCREASING INTENSITY FASTER

If you wish to increase the intensity faster, press the two arrows (▲ fig. C) and the information (i) button simultaneously. This will increase the intensity for both sides at a faster rate.

### Step 4: SESSION OVER

At the end of the training session your THE FLEX BELT by BMR will stop automatically. However, to switch off your unit at any time during a session, press and hold the on/off button for two seconds. You should see the display turn off.

**Step 5:** Remove the belt and lay it on a flat surface. Put the covers back on the black side of the pads. Store your THE FLEX BELT by BMR safely in the pouch until your next session.

**NOTE:** Do not disconnect the unit from the belt without first switching the unit off, as this may result in an error on your display.

**Congratulations on completing your first training session!**

**NOTE:** We recommend that you do at least five training sessions per week to see the best results. You should, however, only do one session per day on a single body area.

## GETTING RESULTS WITH THE FLEX BELT BY BMR

Section 5 provides guidelines for getting the best results as well as some detailed unit information.

1. Use your THE FLEX BELT by BMR 5 times a week.
2. Use the strongest training intensity you can - the harder you work your muscles, the better the results you will see. But remember, stimulation should never be uncomfortable.
3. Use your THE FLEX BELT by BMR as part of a normal, healthy lifestyle. Improve your diet and try to increase the amount of exercise you do.
4. THE FLEX BELT has 10 toning programs which are listed in the table below. There are 8 core programs and 2 bonus programs. Program 1 is set when you switch your unit on for the first time. The unit will then automatically progress through each program up to program 6. Programs 7 and 8 must be selected manually using the program button. The display will show the program currently selected.

THE TABLE BELOW GIVES YOU INFORMATION ABOUT THE 10 TONING PROGRAMS.

CORE PROGRAMS				
Program Name	Session Duration	Training Level	Toning Power	Auto Progression
Beginner	20:00	light	14 $\mu$ C	2 sessions
Intermediate	25:00	moderate	14 $\mu$ C	3 sessions
Advanced	30:00	moderate	16 $\mu$ C	5 sessions
Expert	30:00	Strong	18 $\mu$ C	10 sessions
Pro (i)	30:00	Strong	19 $\mu$ C	10 sessions
Pro (ii)	30:00	Strong	21 $\mu$ C	Indefinite
Ab-Power	30:00	Strong	22 $\mu$ C	Indefinite
Endurance	40:00	Light	14 $\mu$ C	Indefinite
BONUS PROGRAMS				
Program Name	Session Duration	Session Duration	Session Duration	Session Duration
Crunch (i)	3 x 10 crunches	Moderate	16 $\mu$ C	Simultaneous
Crunch (ii)	3 x 10 crunches	Moderate	16 $\mu$ C	Alternate channels

NOTE: All micro-coulomb ( $\mu$ C) values are rounded off to the nearest integer value. This measures the effectiveness of the muscle workout - the higher the  $\mu$ C value, the more effective the muscle workout.

All programs have "warm up" and "warm down" phases. The warm-up prepares your muscles for the main part of your session. The warm down relaxes your muscles as the session is finishing - as you would when exercising.

NOTE: You cannot change a program during a session. You must first switch your unit off and then on again. Then you can select a different program by pressing the program button. Program 7 delivers a strong abdominal workout, which is useful for those involved in sports. Program 8 is lower in intensity and is useful when you require a less vigorous abdominal muscle workout.

### 5. Use THE FLEX BELT by BMR wherever it suits you.

You can use THE FLEX BELT by BMR almost any time and any place. You can use it at home watching television, working at a desk, in the garden or even walking outdoors. THE FLEX BELT by BMR really is a convenient way of toning your abs.

6. To help you get the most from your THE FLEX BELT FROM BMR, we recommend that you follow our 30 day plan. To follow this plan, use your belt 5 times a week for 4 weeks and record your highest intensity levels in your diary (pages 18-19). Each week you have two rest days to allow your muscles to recover. The diary below shows one person's training intensity, recorded during scientific testing over a 30-day period.

NOTE: The intensity level will be different from person to person.

We recommend that you follow a similar plan, trying to push yourself to use higher training intensity levels all the time. Try to increase your training intensity from week to week. But remember, do not over-exert yourself when using muscle stimulation. Any workout should be at a level comfortable for you.

#### 30-Day Plan

Week	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
1	(40 40)	rest	(40 40)	rest	(70 70)	(70 70)	(70 70)
2	rest	(70 70)	(80 80)	rest	(80 80)	(80 80)	(80 80)
3	rest	(85 85)	(85 85)	rest	(85 85)	(90 90)	(150 150)
4	rest	(90 90)	(90 90)	rest	(90 90)	(90 90)	(150 150)

Each box shows the average training intensity from either side of the belt.

Fill-in your diary (pages 18-19) every time you complete a training session. Leave this diary somewhere you will see it regularly. This will help to motivate you to reach your goal. You will soon see and feel that your ab muscles are firmer and more toned.

### CARING FOR YOUR BELT

Your belt can be washed, but you must first remove the unit and pads. Always follow the instructions on the label when washing the belt.

Never machine wash your belt. Always hand wash in lukewarm water. This protects the internal wiring from damage in the washing machine. You should take care even when hand washing and never wring the belt to remove water.

Do not use bleach when washing the belt.

Do not dry clean your belt.

Do not tumble dry your belt. Dry the belt on a flat surface. Do not dry it over anything hot. (e.g. a radiator) as the belt contains plastic parts. Ensure the belt is completely dry before using it again.

The belt should not be ironed.

#### Belt Materials:

Outer material: 100% Nylon, Binding: 82% Nylon, 18% Elastane, Hook & Loop: 100% Nylon, Foam: 100% Polyurethane

### Information Button (i)

Press the information button any time during a session to see the highest intensity level you have reached for that session. Pressing the information button twice shows you the total number of sessions you have completed (Fig. a). Pressing the information button three times shows the belt type you are using.

Fig. a

### Mute Function (M)

If you want to switch off your unit's sound effects, press and hold the program button for two seconds (Fig. b). The mute function remains active indefinitely unless manually changed. Deactivate the mute function by again pressing the program button for two seconds.

Fig. b

### Keylock Function (L)

If you find a comfortable exercise intensity, press and hold the information button for two seconds to lock that intensity (Fig. c). This function only remains active during the session in which it is activated. If you wish to increase the intensity further, deactivate the keylock function by again pressing the information button for two seconds.

Fig. c

### Error Messages

In the unlikely event of your unit developing a problem, "Err" will appear on your unit display (Fig. d). If this occurs you should switch the unit off and then switch it back on again. It should now operate properly. If the problem persists, please call the THE FLEX BELT Careline for further assistance:

Fig. d

**(855) FLEXGEAR (855) 353 - 9432**

E-mail: [info@flexgear.com](mailto:info@flexgear.com)

### Battery Power / Replacing the battery

The battery icon on the display will flash when the battery is low and needs to be recharged. After a period of time, you may find that your battery will only power your unit for a limited number of sessions, requiring you to charge the battery more frequently. This is because rechargeable batteries have a limited number of charge cycles and may eventually need to be replaced.

Fig. e

Should you need to replace the rechargeable battery, remove the small rubber plug from the back cover of your unit (Fig. e), unscrew the back cover and remove it. Replace the existing battery pack with a new battery pack (Fig. f) and replace the battery cover. You can purchase a new battery pack by contacting the THE FLEX BELT by BMR Careline.

Fig. f

## **ADDITIONAL UNIT FUNCTIONS**

### **Battery Charging and discharging**

Your device is powered by a rechargeable battery. The full performance of a new battery is achieved only after two or three complete charge and discharge cycles. The battery can be charged and discharged hundreds of times, but it will eventually wear out. To ensure optimum battery performance you should use your device regularly e.g. perform 1 session a day, 5 days a week.

If a replacement battery is being used for the first time or if the battery has not been used for a prolonged period, it may be necessary to connect the charger and then disconnect and reconnect it to start the charging.

Unplug the charger from the electrical plug and the device when not in use. Do not leave fully charged battery connected to a charger, since overcharging may shorten its lifetime. If left unused, a fully charged battery will lose its charge over time.

If the battery is completely discharged, it may take a few minutes before the charging indicator appears on the display. Use the battery only for its intended purpose. Never use any charger or battery that is damaged.

Leaving the battery in hot or cold places, such as in a closed car in summer or winter conditions, will reduce the capacity and lifetime of the battery. Always try to keep the battery between 15°C and 25°C (59°F and 77°F). A device with a hot or cold battery may not work temporarily, even when the battery is fully charged. Battery performance is particularly limited in temperatures well below freezing. Do not dispose of batteries in a fire as they may explode. Batteries may also explode if damaged. Dispose of batteries according to local regulations. Please recycle when possible. Do not dispose as household waste.

Section 6 has two parts. Part 1 answers frequently asked questions. Part 2 is a problem solving guide.

## HELP PART 1 - FREQUENTLY ASKED QUESTIONS

### **Can I use my unit for post natal exercises?**

- Yes, but you must wait a minimum of **6 weeks** after childbirth before you begin using it and you **must** consult your doctor first.
- If you've had a Caesarean in the past **3 months**, consult your doctor for approval before using the belt.

### **The stimulation is uncomfortable. How can I improve this?**

- Make sure the pads are positioned correctly and that they are pressed firmly against your skin - see page 7 & 8. Switch off your unit and reposition the belt if necessary.
- You can also spread a few drops of water on the black surface of the pads. This can improve the comfort of the muscle stimulation, but be careful not to get water on the unit. Ensure the unit is switched OFF before you do this!
- Ensure the metal studs are fully covered by the pads.

### **I can feel a tingling sensation in my legs during an exercise. What should I do?**

This indicates that the pads are over your hip-bones. Moving the two smaller pads upwards and inwards on your waist (i.e. towards the center of your body) should prevent this. Remember to pause or switch off your THE FLEX BELT by BMR before adjusting the belt or pads.

### **My skin is red after the exercise. Is this a problem?**

- Some redness of the skin after a training session is normal. It is partly due to an increase in the blood flow under the skin and should fade after a while. You may also experience some reddening of the skin due to the pressure of the belt. This is the same as the pressure marks you can get from tight clothing. You should not be concerned about this. It should fade soon after you remove the belt.
- If the redness is excessive, you may have the training intensity too high. This may increase the reddening in sensitive skin. Try using a lower training intensity for a few days. If the problem persists, you should stop using the unit.

### **Will the product cause muscle soreness?**

As with all exercise some muscle soreness can occur after using the belt. This is normal and should go away after a day or two. Use a low training intensity for a few sessions if you are experiencing muscle soreness.

### **How do I know when to replace the pads?**

- With time the pads pick up skin debris and may need to be replaced as this makes the workout less effective and less comfortable.
- You may notice the signal is weakening even if the batteries are OK. This usually indicates that the pads are wearing and will soon need replacing. You can check this by first adjusting the belt to ensure correct positioning of the pads. If the problem persists, try using the unit with a new battery pack. If the signal is still weak or uncomfortable, you should order new pads.
- New pads can be purchased online at [www.theflexbelt.com](http://www.theflexbelt.com) or from the THE FLEX BELT Careline.

### **I can feel my waist muscles exercising but not my stomach muscles.**

Pause the program and reposition the central pad slightly lower on your stomach. If this doesn't help, move the 2 small pads to a smaller figure setting on the belt (towards the central pad).

### **I can feel my stomach muscles exercising but not my waist muscles.**

Pause the program and reposition the central pads slightly higher on your stomach. If this doesn't help, move the 2 smaller pads to a larger figure setting on the belt (away from the central pad).

### **Battery performance has degraded significantly.**

After an extended period of time you may notice a degradation in the performance of the battery in your unit. At this point you should purchase a new rechargeable battery pack. New rechargeable battery packs can be purchased through the THE FLEX BELT Careline, or by visiting [www.theflexbelt.com](http://www.theflexbelt.com).

Problem	Possible Cause	Solution
The display doesn't come on and there is no signal from the unit	The battery is low	Recharge the battery
The display is on but there is no signal	The battery is low	Recharge the battery
The battery symbol is flashing	The battery is low	Recharge the battery
The ⚠ symbol has appeared on the display*	The unit is not properly connected to the belt	Ensure the unit and belt are properly connected
	The covers are still on the pads	Remove the covers from the pads
	The pads are worn	Visit <a href="http://www.theflexbelt.com">www.theflexbelt.com</a> or contact Customer Careline for new pads
	The pads are not covering the studs	Make sure the pads fully cover the metal studs
Unpleasant feeling beneath the pads	Poor pad contact with the skin	Ensure the pads are pressed firmly against the skin. Spread a few drops of water on the surface of the pads. Make sure the unit is OFF first
	The pads are worn	Visit <a href="http://www.theflexbelt.com">www.theflexbelt.com</a> or contact Customer Careline for new pads
	Pads are not covering the metal studs	Reposition the pads so that the metal studs are covered
	Too many consecutive sessions	You should just do one session per day on any single body area
The contractions are very weak even when the training intensity is high	The pads are worn	See <a href="http://www.theflexbelt.com">www.theflexbelt.com</a> or contact Customer Careline for new pads
	The battery is low	Recharge the battery
	Poor pad positioning	See page 7 & 8 for correct pad positioning
	Pads are not covering the studs	Reposition the pads so that the metal studs are covered
"Err" appears on the display	Error message	Switch the unit on and off. The unit should work without any problem. If the error message persists call the Customer Careline for further details

\* The ⚠ symbol always appears in conjunction with the left/ right intensity indicators. If the left indicator appears, there is a problem with the left pad, the right indicator acknowledges a problem with the right pad and if both indicators appear there is a problem with the either the middle pad or all three pads.

## HOW DOES MUSCLE STIMULATION TECHNOLOGY WORK?

Section 7 explains the science behind muscle toning technology.

### How does muscle stimulation technology work?

THE FLEX BELT by BMR uses electro muscle stimulation (EMS) technology. This technology has been used in hospitals and by physical therapists for over 50 years. Here is a summary of how EMS technology works. The special THE FLEX BELT by BMR 3-pad design ensures that all of the abdominal muscles contract and not only those directly under the pads.

**Step 1**  
Signals are sent between the pads

**Step 2**  
The signals switch on the nerves which control the abdominal muscles.

**Step 3**  
The nerves, in turn, cause your muscles to contract and relax rhythmically.

THE FLEX BELT by BMR has all the components of an effective abdominal toner.

1. Your unit has  $21\mu\text{C}$  (micro coulombs) of toning power. This allows you to generate strong, effective muscle contractions.
2. THE FLEX BELT by BMR uses large, medical-grade pads. These ensure that your workout is comfortable.
3. The toning belt's EMS technology has been clinically tested and is proven to give you firmer, stronger, more toned abdominal muscles in 8 weeks.
4. Your unit has 10 programs which cater to a whole range of training levels from the beginner to advanced users.

**NOTE:** Your toning belt will not help you lose weight. For this you need to do some form of cardiovascular exercise e.g. walking, cycling, swimming, jogging. However using your toning belt regularly will give you stronger, firmer, and more toned abs.

## TECHNICAL SPECIFICATIONS

Section 8 is where you will find advanced technical information on THE FLEX BELT BY BMR.

### Caring for your unit

Your unit should not be allowed to get wet or be left in excessive sunlight. It may be cleaned regularly using a soft cloth, lightly dampened in soapy water. Do not allow the interior of your unit to become wet. Do not use detergents, alcohol, spray aerosols or strong solvents on your unit.

Access to the interior of the unit is not required for maintenance purposes.

If your unit is damaged, you should not use it but should return it to THE FLEX BELT or your local distributor for replacement or repair. Repairs, service and modifications may not be carried out by anyone other than qualified service personnel authorized by THE FLEX BELT.

**IMPORTANT:** Under no circumstance should anything other than the correct type of batteries – (rechargeable batteries 3.6V NiMH) be used with your unit. These can be purchased from the THE FLEX BELT Careline.

### Disposing of the pads and batteries

Used pads and batteries must never be disposed of in a fire but in accordance with your country's national laws governing the disposal of such items.

**NOTES:** The signal gradually increases to a peak intensity level at the start of the contraction phase (i.e. ramp up) and gradually decreases to nothing at the start of the relaxation phase (i.e. ramp down). When the  symbol is flashing, the stimulus is significantly reduced.

### Accessories:

You can purchase all accessories at [theflexbelt.com/system](http://theflexbelt.com/system). Under no circumstances should anything other than THE FLEX BELT accessories be used with THE FLEX BELT by BMR charger or THE FLEX BELT by BMR unit (Type 390). Any others may not be compatible with your unit and could degrade the minimum safety levels.

THE FLEX BELT Female Garment (X:70)

THE FLEX BELT Male Garment (X:70)

THE FLEX BELT Pads: • 1 large pad (Type 706 or 709)  
• 2 small pads (Type 707 or 710)

THE FLEX ARMS Garments (E:60)

THE FLEX ARMS Pads: • 4 small pads (Type 715)

THE FLEX MINI Garments (E:30)

THE FLEX MINI Pads: • 4 large pads (Type 716)

**Battery pack** (3.6V, NiMH)

**Battery Charger** (EU - 2504-0303)

**Battery Charger** (US/Japan - 2504-0302)

Charger complies to IEC 60950 and UL 1950, Input: nominal voltage 100-240V, frequency 50-60Hz, current 90mA. Output: 6V / 550mA. Do not use any other chargers or any other power supplies.

**Product Type:** 390

**Intended use:** Muscle stimulator

**Classification:** Internally powered equipment, Class II charger, Type BF applied parts.

**Waveform:** Symmetrical bi-phasic square waveform when measured into a resistive load.

### Environmental Specifications:

**Operating:** Temperature Range: 0 - 35°C (32 - 95°F)  
Humidity: 20 - 65 % RH

**Storage:** Temperature Range: 0 - 35°C (32 - 95°F)  
Humidity: 20 - 85 % RH

### Description of your unit's symbols:

There are a number of technical markings on your unit. These can be explained as follows:

The unit and belt are manufactured for Bio-Medical Research Ltd, Parkmore Business Park West, Galway, Ireland.

The unit requires 1 x 3.6 volt (NiMH) DC battery pack DC is indicated by the symbol: 

The Output Frequency indicates the number of pulses per second transmitted by the unit. This is measured in hertz, which is indicated by "Hz".

Output (RMSA) means the maximum output root mean square current for each channel.

Output (RMSV) means the maximum output root mean square voltage for each channel.

This symbol means "Attention, consult the accompanying documents".

This symbol means type BF applied parts.

SN stands for 'serial number'. Inside the box, on the back of the unit is the serial number specific to this unit. The letter preceding the number indicates the year of manufacture, where "L" denotes 2006, "M" denotes 2007 etc.

The belt's batch number is represented on the belt packaging, by the number corresponding with the LOT symbol.

### Rated Outputs – Voltage/Currents: E :10 and X :10

Parameter	500Ω	1KΩ	1K5Ω
Output RMSV	7.35V	13.1V	15.5V
Output RMSA	14.7mA	13.1mA	10.3mA
Output Frequency	50-70Hz	50-70Hz	50-70Hz
DC Component: Approx.	0 C	0 C	0 C
Positive Pulse Width:	200-315µs	200-315µs	200-315µs
Negative Pulse Width:	200-315µs	200-315µs	200-315µs
Interphase Interval:	100 µs	100 µs	100 µs

### Attention:

This unit can deliver currents in excess of 2mA/cm<sup>2</sup> if used with incorrect electrodes.

At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Some product materials can be re-used if you bring them to a recycling point. By re-using some parts or raw materials from used products you make an important contribution to the protection of the environment. Please contact your local authorities if you need more information about collection points in your area.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

## PRODUCT WARRANTY

Section 9 is where you will find information about our one year limited warranty.

Dear Customer,

Thank you for purchasing your BMR product ("Product"). To be covered by this Limited Warranty, the Product must have been purchased directly from BMR Distribution, Inc. ("SDI"), or its authorized retailer in the United States. This Limited Warranty extends only to the original purchaser and cannot be transferred to anyone. Please read this carefully before using the Product.

### Our Limited Warranty to you

SDI WARRANTS TO YOU, THE ORIGINAL PURCHASER, THAT UNDER NORMAL USE AND SERVICE, THE PRODUCT WILL BE FREE FROM DEFECTS IN MATERIAL AND WORKMANSHIP FOR A PERIOD OF TWO (2) YEAR FROM THE DATE OF PURCHASE ("WARRANTY PERIOD"). IF YOU DISCOVER A DEFECT COVERED BY THIS LIMITED WARRANTY DURING THE WARRANTY PERIOD, AND RETURN THE PRODUCT TO US, WE WILL REPAIR OR REPLACE THE PRODUCT, AT OUR OPTION, USING NEW OR REMANUFACTURED COMPONENTS. THERE IS NO WARRANTY OR CONDITION OF ANY KIND AS TO ANY DEFECTS DISCOVERED AFTER THE WARRANTY PERIOD EXPIRES.

### Certain Product failures not covered

This Limited Warranty covers defects in manufacture that arise from correct use of the Product per the product manual and instructions. It does not cover damage caused by abuse, misuse, accidents, modification or repair by any person other than SDI, moisture, extreme heat or cold, corrosive environments, high voltage surges from external sources or an unsuitable power supply. This Limited Warranty does not cover normal wear and tear on consumable items (e.g. pads, garments, etc.). This Limited Warranty does not apply to any battery provided with the Product, or if the Product has an altered or defaced serial number.

### How to obtain service under this Limited Warranty

To obtain service under this Limited Warranty, please contact SDI's Customer CareLine at 1-800-551-2443 to obtain return instructions. You must submit proof of purchase in the form of a bona fide, dated bill of sale, receipt, or invoice (or copy thereof) evidencing that your request for service is made within the Warranty Period.

Except as may otherwise be required by legislation in your State, SDI will not be responsible for the cost of returning the Product or for any damage to or loss of the Product while it is in transit to us. We recommend that you use a carrier that will track your shipment and ensure the Product's safe delivery. SDI cannot accept any returns sent C.O.D.

**THIS LIMITED WARRANTY IS THE ONLY WARRANTY FOR THE PRODUCT, AND THERE ARE NO OTHER EXPRESS WARRANTIES, ORAL OR WRITTEN, PROVIDED BY SDI. ANY IMPLIED WARRANTIES APPLICABLE TO THE PRODUCT, INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PUR-**

**POSE, ARE LIMITED TO THE DURATION OF THE WARRANTY PERIOD SET FORTH ABOVE.**

### Limitation of Remedies

**YOUR EXCLUSIVE REMEDY, IN LIEU OF ALL INDIRECT, INCIDENTAL, CONSEQUENTIAL, EXEMPLARY, SPECIAL OR OTHER DAMAGES OF ANY KIND, INCLUDING FOR NEGLIGENCE, IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE PRODUCT, IF DEFECTIVE, UNDER THE TERMS AND CONDITIONS OF THIS LIMITED WARRANTY.**

SDI will bear no other expenses.

SDI reserves the right to refuse to repair or replace the Product if it is returned in an unsanitary condition. SDI reserves the right to substitute a product of equal or better quality if an identical product is not available at the time of replacement. This Limited Warranty is non-transferable.

### Applicability of state law

Some states do not allow limitations on the duration of implied warranties, or the exclusion or limitation of incidental or consequential damages, and so the related limitations and exclusions stated above may not apply to you.

### Governing law

This Limited Warranty is valid in the United States only and is not applicable in other countries. If you acquired the Product in the United States, the laws of the State of New York will apply to this Limited Warranty.

**THIS LIMITED WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE.**

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Copyright: © 2012, Bio-Medical Research Ltd. All Rights Reserved
Bio-Medical Research Ltd., Parkmore Business Park West, Galway, Ireland.