Sunbeam

Back Contouring Heating Pad

Model: 300

IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

DANGER: TO REDUCE RISKS OF BURNS, ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

1. READ ALL INSTRUCTIONS.
2. DO NOT USE WHILE SLEEPING.
3. DO NOT USE ON AN INFANT OR ON AN ANIMAL.
4. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, HELPLESS, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION.
5. DO NOT USE ON AREAS OF INSENSITIVE SKIN.
6. BURNS CAN OCCUR REGARDLESS OF HEAT CONTROL SETTING. CHECK SKIN UNDER PAD FREQUENTLY.
7. NEVER USE PAD WITHOUT COVER IN PLACE.
8. DO NOT USE IN AN OXYGEN ATMOSPHERE.
9. DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.
10. DO NOT BUNCH, SIT ON OR CRUSH PAD-AVOID SHARP FOLDS.
11. NEVER PULL THIS PAD BY THE SUPPLY CORD.
12. DO NOT USE THE CORD AS A HANDLE.
13. CAREFULLY EXAMINE INNER COVER BEFORE USING. DISCARD THE PAD IF INNER COVER SHOWS SIGN OF DETERIORATION.
14. THIS UNIT SHOULD NOT BE USED BY CHILDREN WITHOUT ADULT SUPERVISION.
15. USE THIS PAD ONLY ON A 110-120 VOLT AC CIRCUIT. UNPLUG WHEN NOT IN USE.
16. DO NOT TAMPER WITH THIS PAD IN ANY WAY.
17. THERE ARE NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PAD DOES NOT FUNCTION SATISFACTORILY, SEE WARRANTY FOR CONSUMER SERVICE INFORMATION.
18. DO NOT USE THIS PAD WITH LINIMENT, SALVE OR OINTMENT PREPARATIONS, ESPECIALLY ONES THAT CONTAIN HEAT-PRODUCING INGREDIENTS. SKIN BURNS COULD RESULT.
19. DO NOT BEND OR PINCH CORD AND DO NOT USE IF CORD IS DAMAGED.
20. NEVER LEAVE THE APPLIANCE UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.
21. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.
22. WHEN NOT IN USE, STORE THE HEATING PAD IN THE ORIGINAL PACKAGING, IN A DARK, COOL, DRY PLACE.
23. SAVE THESE INSTRUCTIONS.

THIS APPLIANCE HAS A POLARIZED PLUG (ONE BLADE IS WIDER THAN THE OTHER). TO REDUCE THE RISK OF ELECTRIC SHOCK, THIS PLUG IS INTENDED TO FIT IN A POLARIZED OUTLET ONLY ONE WAY. IF THE PLUG DOES NOT FIT FULLY IN THE OUTLET, REVERSE THE PLUG. IF IT STILL DOES NOT FIT, CONTACT A QUALIFIED ELECTRICIAN TO INSTALL THE PROPER OUTLET.

DO NOT ATTEMPT TO MODIFY THIS PLUG OR DEFEAT THIS SAFETY FEATURE IN ANY WAY.

HOW TO USE YOUR BACK CONTOURING HEATING PAD

1. In addition to the therapeutic heat, this pad features a self-inflatable lumbar support, which can be adjusted for designed firmness. To INFLATE the pad, open the Air Inlet Valve (round black knob on side of pad) and the pad will self inflate automatically. You can regulate the firmness of the pad by closing the air valve when the pad reaches the desired thickness. To DEFLATE the pad, open the air inlet valve and roll up the pad, forcing the air out. The pad can be stored in this rolled up position.
2. Plug into any 120V AC electrical outlet.
3. Place the pad on a chair or other similar location where you plan to use it.

NOTE: For optimum heat therapy, ensure that the pad's softer side with the "Contouring Airflow System" logo is positioned towards the body. The pad is designed to have the heated surface only on one side of the pad. Follow the markings on pad's cover for more direction.

WARNING: Never use this pad without the cloth cover.

4. To turn the pad on, press the Power Button and set control to the desired heat setting. To turn the pad OFF, scroll through remaining heat settings by pressing the Power Button until none of the heat settings are selected.
5. Auto-Off function will shut the heating pad down automatically after 30 minutes of use. This is indicated by the flashing light of Auto-Off indicator on the controller. To reset, press the Power Button and select the desired heat setting.

NOTE: In case of power outage, the Auto-Off feature will automatically reset the timer.

USING MOIST HEAT THERAPY

In some cases, you may find that moist heat therapy provides more comfort and relief. **When using moist heat therapy, you should protect furniture and clothing from becoming damp.**

1. Dampen one side of the cover with warm water by patting the surface with a wet sponge or cloth or use a spray bottle to apply a fine mist. Do not thoroughly soak the cover.
2. Do not spray or dampen the power cord.

DANGER: Do not immerse the pad or wiring in water or place under running water. Failure to comply may result in product failure and/or personal injury.

FABRIC COVER WASHING INSTRUCTIONS

DO NOT DRY CLEAN! DO NOT BLEACH! DO NOT USE WRINGER!

1. Disconnect pad from outlet.
2. Remove cloth cover from pad.
3. Machine-wash the cover in cool water on gentle cycle, or hand wash. Use mild detergent.
4. Hang to dry or machine dry on gentle, cool cycle.
5. Insert the vinyl pad inside the cloth cover.

NOTE: The cover is designed so that the air inlet valve fits inside it only one way. This feature ensures that the pad's heated surface is always on the soft side of the cover, which is intended to be used against the body. Follow the markings on the cover for more direction. Make sure that the cord and controller remain on the outside of the cover. Zip the cover closed.

WARNING: Never use this pad without the cloth cover.

5-YEAR LIMITED WARRANTY

Sunbeam Products, Inc. doing business as Jarden Consumer Solutions or if in Canada, Sunbeam Corporation (Canada) Limited doing business as Jarden Consumer Solutions (collectively "JCS") warrants that for a period of five years from the date of purchase, this product will be free from defects in material and workmanship. JCS, at its option, will repair or replace this product or any component of the product found to be defective during the warranty period. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. This is your exclusive warranty. Do NOT attempt to repair or adjust any electrical or mechanical functions on this product. Doing so will void this warranty. This warranty is valid for the original retail purchaser from the date of initial retail purchase and is not transferable. Keep the original sales receipt. Proof of purchase is required to obtain warranty performance. JCS dealers, service centers, or retail stores selling JCS products do not have the right to alter, modify or any way change the terms and conditions of this warranty. This warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, disassembly, repair or alteration by anyone other than JCS or an authorized JCS service center. Further, the warranty does not cover: Acts of God, such as fire, flood, hurricanes and tornadoes.

What are the limits on JCS's Liability?

JCS shall not be liable for any incidental or consequential damages caused by the breach of any express, implied or statutory warranty or condition. Except to the extent prohibited by applicable law, any implied warranty or condition of merchantability or fitness for a particular purpose is limited in duration to the duration of the above warranty. JCS disclaims all other warranties, conditions or representations, express, implied, statutory or otherwise. JCS shall not be liable for any damages of any kind resulting from the purchase, use or misuse of, or inability to use the product including incidental, special, consequential or similar damages or loss of profits, or for any breach of contract, fundamental or otherwise, or for any claim brought against purchaser by any other party. Some provinces, states or jurisdictions do not allow the exclusion or limitation of incidental or consequential damages or limitations on how long an implied warranty lasts, so the above limitations or exclusion may not apply to you. This warranty gives you specific legal rights, and you may also have other rights that vary from province to province, state to state or jurisdiction to jurisdiction.

How to Obtain Warranty Service

In the U.S.A. - If you have any question regarding this warranty or would like to obtain warranty service, please call 1-800-435-1250 and a convenient service center address will be provided to you.

In Canada, - If you have any question regarding this warranty or would like to obtain warranty service, please call 1-800-667-9623 and a convenient service center address will be provided to you.

In the U.S.A., this warranty is offered by Sunbeam Products, Inc. doing business as Jarden Consumer Solutions located in Boca Raton, Florida 33431. In Canada, this warranty is offered by Sunbeam Corporation (Canada) Limited doing business as Jarden Consumer Solutions, located at 20 B Herford Street, Brampton, Ontario L6Y 0M1. If you have any other problem or claim in connection with this product, please write our Consumer Service Department.

PLEASE DO NOT RETURN THIS PRODUCT TO ANY OF THESE ADDRESSES OR TO THE PLACE OF PURCHASE.

©2012 Sunbeam Products, Inc. doing business as Jarden Consumer Solutions. All rights reserved. Distributed by Sunbeam Products, Inc. doing business as Jarden Consumer Solutions, Boca Raton, Florida 33431.
Printed in China. DW-052912

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