

Pop Sonic UV Sanitizing Toothbrush – Tips & Tricks

Sanitizing:

- The Pop Sonic UV Case automatically emits a blue light to sanitizer your toothbrush head when the case closes.
- The UV light will work without being plugged in for the ultimate convenience.
- Any Pop Sonic Toothbrush head can be used with this toothbrush & case

Charging:

- The Sonic Toothbrush charges inside the case for convenience. Put the toothbrush away in its case and plug it in!
- Pop Sonic recommends 16 hours of charging the case with the toothbrush inside before use for maximum sonic power.
- The case and the toothbrush each hold their own charge.
- When you plug the case in, the “charge” light on the toothbrush will light up. When the toothbrush is fully charged, the charge light will turn OFF.

