

NuFACE TRINITY® ELE ATTACHMENT

The dual spheres of the NuFACE Trinity ELE (Trinity ELE) are similar to that of the NuFACE Trinity Facial Trainer (Facial Trainer) in that they are capable of providing complete facial stimulation. However, the dual spheres of the Trinity ELE are significantly smaller in size and designed for optimal contact with smaller surface areas of the face, like the areas around the lips and eyes. The Trinity ELE is extremely responsive, as it delivers low-level microcurrent and allows users to adjust the output level for personalized comfort.

IDEAL USER

The Trinity ELE is for anyone looking to achieve facial stimulation for strategic locations on the face, mainly around the lip and eye areas.

GETTING TO KNOW YOUR DEVICE

1. POWER ON/OFF BUTTON
2. INTENSITY CONTROL UP
3. INTENSITY CONTROL DOWN
4. TRINITY ELE ATTACHMENT
5. SPHERE
6. ATTACHMENT RELEASE BUTTON

