

PMD Replacement Discs



White Disc

Training / Ultra Sensitive



For people with ultra sensitive skin or those first learning to use the PMD.



Blue Disc

Sensitive



For people with sensitive skin that have used the PMD for several treatments and are ready for a level of exfoliation above the Grey Disc.



Green Disc

Medium



For people with normal or non-sensitive skin that have used the PMD for several treatments and are ready for a level of exfoliation above the Sensitive Disc.



Red Disc

Coarse



For people with normal to very coarse skin that have used the PMD for an extended period of time and would like an intense level of exfoliation.



FEEL GREAT IN YOUR OWN SKIN