



HairMax® Professional 12 PARTS


TOP VIEW HOUSING



BOTTOM VIEW HOUSING



SIDE VIEW HOUSING



PARTS:

- (a) BATTERY INDICATOR - Yellow during charging. Green is fully charged.
- (b) INFORMATION DISPLAY - Laser ON, Battery Level, Usage Time.
- (c) POWER BUTTON with GREEN light indicator - Turns on power.
- (d) LASER BUTTON with RED light indicator - Turns on laser (POWER button must be pressed first).
*Caution - Use only as directed. Any variation of controls, adjustments, performance or procedures, other than those specified, may result in hazardous exposure to laser light.
- (e) SAFETY LABEL (removable) - "PLEASE REFER TO OUR USER MANUAL BEFORE USING THIS PRODUCT."
- (f) POWER CORD
- (g) TEETH - removable for cleaning
- (h) LASER WINDOW - Do not look into laser window when laser is turned on. See CAUTIONS
- (i) LASER WARNING LABEL - contains SERIAL NUMBER - not removable.
- (j) CHARGING CRADLE (STAND) - Stores and charges HairMax.
- (k) CHARGING CRADLE INDICATOR - RED during charging, GREEN is fully charged.