

TAKING CARE OF YOUR TRIA AGE-DEFYING LASER

CLEANING

For best results, keep the laser treatment tip clean. Oils, lotions, makeup, or other debris can block the laser light or keep the device from functioning properly.

If the laser treatment tip does become visibly dirty, clean it with a lint-free cloth.

If other parts of the device or charging cradle become dirty, clean as needed with a soft cloth moistened with rubbing alcohol.

Make sure the charging cradle is not plugged in when cleaning. Avoid getting water or other liquid inside any part of the device or the charging cradle.

MAINTENANCE

The device does not need regular or special maintenance. If you think your device or charging cradle is not working properly or is damaged, do not use the device. Call Tria Beauty Customer Care for assistance.

STORAGE

Keep the device and charging cradle away from water, dampness, or debris that may get inside the device or the charging cradle.

DISPOSAL

This device contains lithium batteries. Please follow the regulations where you live for recycling or disposing of electrical devices.

TRAVEL

If you travel with your Tria Age-Defying Laser, pack the device and charging cradle in its original packaging to prevent damage.

The charging cradle will operate properly at all common household voltages (any supply between 100 and 240 VAC at 50-60 Hz). You may need a plug adapter so that the charging cradle can plug into the wall socket.