

product care and maintenance

Please keep the lens / stainless steel retainer assembly clean and avoid scratching the lens as the LED light may be reduced in efficiency which can lead to a decrease in effectiveness of treatment.

SKINCLINICAL reverse's plastic housing may be cleaned ONLY with a clean soft dry NON-ABRASIVE cloth OR an Isopropyl Alcohol wipe. Please follow these guidelines as cleansers or other more caustic materials may damage both the housing and the lens surface.

The wall charger unit should not require any cleaning. If you choose to do so, ALWAYS unplug it and ONLY wipe the outside surfaces with a DRY cloth.