

FAQs

1. How often can I use GloPRO? Do I have to use it at night?

For best results, GloPRO should be used at least 2-3 times a week, however it can be used daily. GloPRO is best used at night but can be used at any time.

2. Which products can I apply to maximize results?

Product absorption is enhanced when applied immediately after GloPRO use. We recommend applying all targeted topicals including serums, oils and moisturizers directly after using GloPRO for best results.

3. How much pressure should I apply with GloPRO?

There is no need to apply additional pressure, just let the weight of the tool work its magic as you glide it over the skin.

4. Why are there so many different MicroTips? Do I need all of them? Each MicroTip Attachment Head is designed to target a specific area of skin. Most GloGetters love them all, but choose what's right for you. The different optimal lengths help address different areas.

5. Is there anywhere on the skin that I should avoid using GloPRO? We do not recommend using GloPRO on the delicate skin on the eyelid. We also recommend avoiding usage of GloPRO on any skin irregularities, which include but are not limited to: acne, open wounds, scars, eyelids or any raised surfaces including moles. If you have any questions or concerns about specific areas, please contact your physician.

6. How often should I replace the MicroTip Attachment heads?

We recommended swapping them out every 3 months, just like you would replace a toothbrush for optimal benefits.

If you have any additional questions, please contact a customer care specialist at 877-945-5274 or info@beautybio.com. For additional information, please visit BeautyBio.com

Information provided by BeautyBio

BEAUTYBIO®