




BEAUTYBIO®

# Getting Started



Information provided by BeautyBio

BEAUTYBIO®

# Care

Follow these steps to clean your GloPRO®:

1. Remove cap from empty spritzer bottle provided.
2. Fill bottle with isopropyl alcohol (70% or higher).
3. Place cap on bottle and tighten.
4. After each use, spray MicroTip™ Head to sterilize the needles and let air dry.



# Storage

When not in use, store your GloPRO® in a clean, dry environment.

**Bonus:** Each replacement GloPRO® MicroTip™ Attachment Head comes with a protective clear cover for easy storage and travel.

Information provided by BeautyBio

BEAUTYBIO®