


When + How



MINIMUM 2-3
TIMES PER WEEK



60 SECONDS PER
TARGET AREA



GLIDE TOOL IN ALL
DIRECTIONS



FACE



EYES

1. Prior to GloPRO use, cleanse skin and prep with GloPRO Skin Prep Pads. The GloPRO Skin Prep Pads can be used on face and neck to help condition the skin before microneedling.
2. Glide GloPRO on target areas applying light pressure rolling in vertical, horizontal and diagonal directions. Spend 60 seconds per target area. Treated areas may appear flushed after use.

#GloTip: Tool can be used in the on or off mode. However, tool must be on to engage LED red light and VibroTactile Technology. Comes with two fully charged AAA batteries. Batteries are replaceable but not rechargeable.

Information provided by BeautyBio

BEAUTYBIO®