

# FAQs

## **1. How often can I use GloPRO? Do I have to use it at night?**

For best results, GloPRO should be used at least 2-3 times a week, however it can be used daily. The GloPRO body roller can be used at any time of day.

## **2. Which products can I apply to maximize results?**

Product absorption is enhanced when applied immediately after GloPRO use. We recommend applying all targeted topicals including body creams, oils and topicals directly after using GloPRO for best results.

## **3. How much pressure should I apply with GloPRO?**

There is no need to apply additional pressure, just let the weight of the tool glide over the skin.

**4. Why are there so many different MicroTips? Do I need all of them?** Each MicroTip Attachment Head is designed to target a specific area of skin. Most GloGetters love them all, but choose what's right for you. The different optimal lengths help address different areas.

**5. Is there anywhere on the skin that I should avoid using GloPRO?** We do not recommend using GloPRO on the delicate skin on and around genitals. We also recommend avoiding usage of GloPRO on any skin irregularities, which include but are not limited to: acne, open wounds, scars, or any raised surfaces including moles. If you have any questions or concerns about specific areas, please contact your physician.

## **6. How often should I replace the MicroTip Attachment heads?**

We recommended swapping them out every 3 months, just like you would replace a toothbrush for optimal benefits.

If you have any additional questions, please contact a customer care specialist at 877-945-5274 or [info@beautybio.com](mailto:info@beautybio.com). For additional information, please visit [BeautyBio.com](http://BeautyBio.com)

Information provided by BeautyBio

BEAUTYBIO®