

Actonis™

USER MANUAL

DUAL ACTION **18**

LaserComb

INTRODUCTION

Introduction to Actonis™

The Actonis Dual Action 18 LaserComb is a home-use laser device that has been clinically proven to treat early to moderate stages of hereditary/hormonal hair loss.

Indication

Actonis models are indicated to treat Androgenetic Alopecia, and promote hair growth in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss with Fitzpatrick Skin Types I to IV.

How Actonis works

Actonis works by providing therapeutic laser energy to stimulate hair follicles. Your device should be used 3 times per week, spaced out every other day. For example, Monday, Wednesday, Friday. The device should not be used more than the recommended treatment time per day. Clinical studies demonstrate that most users see benefits in as little as 16 weeks.

Actonis Results

Actonis laser therapy has been proven in clinical studies. In these clinical studies, over 97% of the women experienced significant hair growth, on average, 21 new hairs per square centimeter/133 per square inch.

It is important that you read and understand the content of this user manual before using your device.


WHO IS THE BEST CANDIDATE FOR Actonis?

How do you know if the Actonis LaserComb is right for you?

Your device can be used by women with thinning hair (including color treated hair) or pattern baldness caused by a hereditary/hormonal condition. Doctors use a system known as the Ludwig-Savin Classification (women) to describe the degree of hair loss. Below are the pictures of the scales. The shaded areas show the type of hair loss that can be treated with your device.

CLASSIFICATION FOR WOMEN

Ludwig-Savin Classification - Females



SYMBOLS

- | | | | |
|--|---------------------------|--|-------------------------------|
| | European Conformity | | Power On/Off |
| | WEEE | | Refer to instruction manual |
| | Warning | | RF interference |
| | Caution | | Type BF Applied Part |
| | Laser | | Serial Number |
| | Laser warning | | Degree of moisture resistance |
| | Laser aperture indication | | Low-level ESD |
| | Laser warning label | | Manufacturer |

TESTING FOR SENSITIVITY

Prior to treatment, test your skin for sensitivity by shining the light on the inside of your forearm for 2 minutes. If you see any reaction, do not use your device and contact customer service.

HOW TO MONITOR DEVICE ACTIVITY

Your device is manufactured in an ISO 13485 certified facility using the highest quality components. Your device should provide many years of trouble-free usage. When working properly, a therapeutic red laser light should shine out of the device. The laser light does not diminish over time.

PARTS



This is a harmless low-level ESD (Electrostatic Discharge) sensitive device. Please do not touch the charging contacts or USB-port. Touching the contacts or USB-port is not of the essential performance and will not cause a basic safety issue. If touching the contacts or USB-port causes the device to malfunction, reset your device by powering off and back on.

WHAT IS INCLUDED

1. Actonis 18 LaserComb
2. Charging cradle
3. Power supply with cord
4. Quick Start Guide

Notice: Entire Actonis assembly is considered as a Type BF applied part.

GETTING STARTED

Your Actonis device comes with a rechargeable battery permanently installed. Read this manual completely before you start to use your Actonis. The instructions in this user manual are intended for individuals with a minimum of an eighth grade education. Actonis is safe for home use. The consumer is the intended operator of this equipment.


Use the Actonis Dual Acton 18 LaserComb for 8 minutes per treatment, three treatments per week. The treatments should be spaced every other day, for example Monday, Wednesday, Friday. Do not use your Actonis more than the recommended treatment time per day. Most users see benefits in as little as 16 weeks (based on clinical studies).




When you receive your device, take all the parts out of the case. You must charge your Actonis for 2-3 hours before initial use.

CHARGING OPTIONS

OPTION 1: Charge your device using the charging cradle and power cord (Figure A). Prior to placing your device in the cradle, make sure it is turned off and the teeth face away from you. Plug the power cord into the cradle and the other end into an electrical wall outlet. A green light will illuminate on the back of the charging stand indicating that your power cord is operating properly.




OPTION 2: Charge your device by placing the power cord directly into the device and the other end into an electrical wall outlet (Figure B).



The indicator light on the power button will flash during the charging process. Once charging is complete, the light will stop flashing and remain on.

When the battery is fully charged, the indicator light will be on. As the battery power diminishes, the light will blink slowly. When the battery is getting low and needs to be recharged, the light will start to flash rapidly. When this happens, charge your device with one of the charging options as previously described.

During use, the indicator light will indicate the battery level. A steady, non-flashing light indicates the battery is full (Figure C).



If your device needs to be recharged during your treatment, allow it to fully charge before continuing.

PROPER USE


The Actonis Dual Action 18 LaserComb should be used for 8 minutes.

Use a watch or clock to time your treatment.

First, press the power button until power button light turns on.



Second, press the laser button. Your device is ready for use.



Shine your device onto the palm of your hand and you will see that the laser light is on (do not look at the laser light directly).


Do not use your device if the device is damaged or defective, contact Actonis customer service or your place of purchase.

If your device needs to be recharged during your treatment, allow the device to fully charge before continuing your treatment. Do not turn on or use your device while it is charging.

TREATMENT METHOD

Prior to treatment, test your skin for sensitivity. See prior TESTING FOR SENSITIVITY section.

After you have turned on your device, place it flat on your scalp so that both rows of the teeth are touching your scalp (Figure D).



Start at your hairline. Move your device from front to back of scalp in a slow combing motion, pausing every 1/2 inch (1 cm) for 4 seconds (Figure E). When your device is on, it will beep and vibrate every four (4) seconds. This means that it's time to move your device to the next area of the scalp. Do not lift your device off your scalp while combing from front to back.

Continue in this manner from front to back and then side to side until you cover your entire scalp with multiple passes (Figure E).



Repeat these passes for the total treatment of 8 minutes.

If you have long hair, use your fingers to lift your hair out of the way while using your device to be sure that the laser energy reaches your scalp.



NOTE: The Actonis Dual Action 18 LaserComb will turn off automatically after 8 minutes.



IMPORTANT: Do not lift your device off your head when moving from front to back. Your device features a pair of comb teeth that part the hair so the light can more effectively reach the scalp. You will make a few complete passes over your scalp within the recommended treatment time.

If you wish to **turn off the beep/vibration** during the treatment, press the top laser button once (Figure G). Your device and laser will remain on. To re-activate beep and vibration, simply press the laser button twice.

Turn off your device by pressing and holding the power button (lower button) until indicator light turns off (Figure H). Store the device in the charging cradle for future use. You may leave the power cord plugged into the charging cradle at all times.



Please Note:

It is important to move your device very slowly. Remember, you are bathing the scalp with laser light. If you move it too fast, the treatment will not achieve the optimal results to promote hair growth.

TROUBLESHOOTING

ANY ATTEMPT TO SELF-SERVICE YOUR DEVICE, UNLESS INSTRUCTED BY ACTONIS CUSTOMER SERVICE, WILL VOID YOUR WARRANTY.

Actonis has advanced electronics to ensure that the device is working like it is supposed to and that the proper amount of laser energy is being delivered. If for some reason you are experiencing difficulties with your device, please try one of the tips below:

If your device will not turn on

- Make sure your device is fully charged.
- Press the power button for 4 seconds and then the laser button.
- Try another electrical wall outlet.

If your device is turning off during use

- Make sure your device is fully charged.
- Press the power button for 4 seconds and then the laser button.
- Try another electrical wall outlet.

If your device will not turn off

- Press and hold the power button for 4 seconds.

If your device beeps/vibrates irregularly or doesn't beep/vibrate at all

- Press the power button then the laser button.

If your device beeps and vibrates, but you do not see the laser light

- Press the power button then the laser button.

LASER PARAMETERS

Wavelength: 655nm and 635nm visible red light

Laser Power for classification: <5mW

Divergence: 57 mrad

Number of lasers: 9 (655nm) 9 (635nm)

This Laser Product is designated as Class 3R during all procedures of operation.

Note: Actonis laser devices have no user-serviceable parts.

USE CONDITIONS

MTBF Lasers (mean time between failures): 10,000 hours at 25°C ambient temperature

Operating temperature/altitude/atmospheric pressure/humidity: 5°C to 35°C/0 to 2,000 m/80kPa to 106kPa/0% to 90%, relative humidity

Storage & transportation temperature/humidity/atmospheric pressure: -20°C to 60° C/0% to 90%, relative humidity/50kPa to 106kPa

Ingress of water or particulate matter: IP22

Battery charging cycles: 300 cycles

ROHS (restriction of hazardous substances): 2011/65/EU

Device size and weight: 9.5 x 2.5 x 2 in (24.3 x 6 x 5 cm) / 6.4 oz (181 gm)

Cradle size and weight: 5.5 x 4.375 x 2.75 in (14 x 11 x 7 cm) / 7.0 oz (198 gm)

Product warranty: 2 years

Age of use: Adult

Software Version: KM910DW-V19-20230613A-120hz

WARNING LABEL

Warning label



Serial number

**POWER SUPPLY SPECIFICATIONS**

The power supply for your Actonis can be used in the United States and Internationally. This power supply will adapt to all power supplies from 100 volt/60 Hz (US) to 240 volt/50 Hz. If outlets in your country have different plugs than the power supply that comes with the Actonis, use an inexpensive plug adapter not supplied with the device.

WARNINGS

- ▲ Do not stare directly at the red laser light or its reflection in a mirror because it can temporarily irritate your eyes. If you stare at the light for too long, it could harm your eyes. Also, never use a magnifying glass when using your device because it can cause temporary irritation to your eyes.
- ▲ Keep the device out of the reach of children and pets.
- ▲ Do not drop your device in water because you may get an electrical shock. If your device gets damaged by water, contact customer service.
- ▲ If the cord is damaged on your device do not use the device because you could get an electrical shock. If the cord gets damaged, contact customer service.
- ▲ Do not use your device near any heated surfaces because this could cause the device not to work properly and may cause an electrical shock, contact customer service.
- ▲ Laser light - avoid direct eye exposure.
- ▲ This is a harmless low-level ESD (Electrostatic Discharge) sensitive device. Please do not touch the USB port. Touching the contacts or USB port is not of the essential performance and will not cause a basic safety issue. If touching the contacts or USB port causes the device to malfunction, reset your device by powering off and back on.
- ▲ No modification of this equipment is allowed.
- ▲ This equipment has no user-serviceable parts. Do not attempt to modify or repair this equipment.
- ▲ Do not attempt to modify or service while the equipment is in use.
- ▲ If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use. Contact Actonis customer service.

CAUTIONS

- △ If you are taking any medications or products that make you sensitive to light, you should test your device to make sure you do not have a reaction. Some of the things that can cause increased sensitivity to your device could be medicines for colds, allergies, pain medicine and drugs to treat infections. To test your device for sensitivity, turn the device on, hold the device on one spot, shining the light on your forearm. After two minutes, turn the device off and set aside. Look at the spot where you held the laser light over the next five minutes. If you see any reaction such as redness or rash, do not use your device and contact our customer service information on how to return to place of purchase and get a refund.
- △ Only use your device as directed. Do not try to change the way your device operates. Altering controls, adjustments, performance or procedures other than those specified herein may result in hazardous radiation (laser) exposure. Contact customer service.
- △ Never use acetone or other solvents on any part of your device. Use of solvents on the window can block the laser light from coming through.

MAINTENANCE

Actonis is designed to last for many years. It does not require any maintenance other than periodic cleaning.

Your device is intended for indoor use only.

CLEANING YOUR ACTONIS

- Keep your device clean.
- Turn off the device by pressing and holding the power button (lower button) until the light turns off.
- Unplug the power cord from the device and from the electrical wall outlet.
- To clean the teeth more thoroughly, slide the teeth towards the top of your device and gently pull out. Clean with soap and water.
- Clean the laser window with a soft, slightly moistened cloth. Don't use any acetone or solvents.

DISPOSAL INSTRUCTIONS



This symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

ELECTRICAL SAFETY

Basic safety precautions should be observed when using this, or any other electrical device, especially when children are present.

▲ WARNING: Keep away from water. Contact with water may result in an electrical shock.

To reduce the risk of hazard by electric shock:

- ▲ Do not use while bathing.
- ▲ Do not place or store your Actonis where it can fall or be pulled into a tub or sink.
- ▲ Do not expose your Actonis to running water, or drop into water or any other liquid. If your Actonis falls into water, unplug it immediately from the power outlet.
- ▲ If using your Actonis near water, the use of an electrical outlet protected by a Ground Fault Circuit Interrupter is recommended.
- ▲ Do not dispose of your Actonis in a fire or store in high temperatures greater than 140°F (60°C).
- ▲ Do not turn on or use your Actonis while it is charging.
- ▲ Power supply cannot be separated during use. Do not position the equipment in a manner that will make it difficult to plug in or remove the power cord.
- ▲ If you have any problems with this device, such as setting up, maintaining or using, please contact customer service.
- ▲ Please report to customer service if any unexpected operation or events occur.

- ▲ After drop/shock, that may cause changes in the performance, please contact customer service.
- ▲ Be careful, do not put the power cord around your neck, it may cause strangulation.
- ▲ Keep the device out of the reach of children/pets to avoid inhalation or swallowing of small parts.
- ▲ Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- ▲ The device should be used only with the accessories recommended for use by the manufacturer.
- ▲ When not in use, store the device with the adapter in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.

ELECTROMAGNETIC EMISSIONS GUIDANCE

Guidance and manufacturer's declaration – electromagnetic emissions

The LaserComb intended for use in the electromagnetic environment specified below. The customer or the user of the LaserComb should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The LaserComb use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The LaserComb is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Compliance	

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) the device must accept any interference received, including interference that may cause undesired operation.

ELECTROMAGNETIC IMMUNITY GUIDANCE

Guidance and manufacturer's declaration – electromagnetic immunity

The LaserComb is intended for use in the electromagnetic environment specified below. The customer or the user of the LaserComb should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2kV for power supply lines ± 1kV for input/output lines	± 2kV for power supply lines Non applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1kV line(s) to line(s) ± 2kV line(s) to earth	± 1kV differential mode Non applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the LaserComb require continued operation during power mains interruptions, it is recommended that the LaserComb be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The LaserComb power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

ELECTROMAGNETIC IMMUNITY GUIDANCE

Guidance and manufacturer's declaration – electromagnetic immunity

The LaserComb is intended for use in the electromagnetic environment specified below. The customer or the user of the LaserComb should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the LaserComb including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the LaserComb is used exceeds the applicable RF compliance level above, the LaserComb should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the LaserComb. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

ELECTROMAGNETIC DISTANCE RECOMMENDATIONS

Recommended separation distance between portable and mobile RF communications equipment and the LaserComb.

The LaserComb is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the LaserComb can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the LaserComb as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

WARRANTY

ANY ATTEMPT TO SELF-SERVICE YOUR DEVICE, UNLESS INSTRUCTED BY LEXINGTON INTERNATIONAL CUSTOMER SERVICE, WILL VOID YOUR WARRANTY

To validate your warranty, please register online at: www.actonis.com/warranty

THE DEVICE IS WARRANTED BY LEXINGTON INTERNATIONAL, LLC ("LEXINGTON") AGAINST DEFECTS IN MATERIALS OR WORKMANSHIP FOR 12 MONTHS FROM DATE OF PURCHASE. IF THIS PRODUCT IS DETERMINED TO BE DEFECTIVE, LEXINGTON WILL REPAIR OR REPLACE IT, AT ITS OPTION, AT NO CHARGE DURING THE WARRANTY PERIOD. THE WARRANTY PERIOD IS NOT EXTENDED IF LEXINGTON REPAIRS OR REPLACES THE PRODUCT. AFTER THE WARRANTY PERIOD, THE USER MUST PAY FOR ALL PARTS AND LABOR CHARGES. IF THE PRODUCT IS REPAIRED OR REPLACED, EITHER WITHIN OR OUTSIDE THE WARRANTY PERIOD, LEXINGTON RESERVES THE RIGHT TO PROVIDE EQUIVALENT PARTS IN DIFFERENT COLORS TO YOUR ORIGINAL PRODUCT.

THIS WARRANTY DOES NOT COVER DAMAGES CAUSED BY ACTS OF GOD, MISUSE, NEGLIGENCE, ACCIDENT, MODIFICATION OF OR TO ANY PART OF THE PRODUCT. THIS WARRANTY DOES NOT COVER DAMAGE DUE TO IMPROPER OPERATION OR MAINTENANCE, CONNECTION TO IMPROPER VOLTAGE SUPPLY, OR ATTEMPTED REPAIR BY ANYONE OTHER THAN LEXINGTON. THIS WARRANTY IS VOID IF PRODUCT IS TAMPERED WITH.

TO OBTAIN WARRANTY SERVICE, YOU MUST CALL LEXINGTON INTERNATIONAL CUSTOMER SERVICE TO REQUEST AN AUTHORIZATION NUMBER. YOU MUST DELIVER THE PRODUCT FREIGHT PREPAID, IN EITHER ITS ORIGINAL PACKAGING OR PACKAGING AFFORDING AN EQUAL DEGREE OF PROTECTION, TOGETHER WITH YOUR NAME, ADDRESS AND AUTHORIZATION NUMBER, TO:

Lexington International, LLC, 1040 Holland Drive, Boca Raton, FL 33487, USA

LEXINGTON INTERNATIONAL, LLC SHALL NOT BE LIABLE FOR ANY SPECIAL, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES THAT ARE DIRECTLY OR INDIRECTLY RELATED TO THE USE OF, OR THE INABILITY TO USE, THE PRODUCT, EVEN IF LEXINGTON INTERNATIONAL, LLC HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES WHETHER IN AN ACTION OF CONTRACT, NEGLIGENCE OR OTHER TORTUOUS ACTION. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF SPECIAL, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU. THIS WARRANTY IS THE ONLY WRITTEN OR EXPRESS WARRANTY GIVEN BY LEXINGTON INTERNATIONAL, LLC. ALL IMPLIED WARRANTIES FOR THE PRODUCT, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN TIME TO THE TERM OF THIS WARRANTY. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU.

GOVERNING LAW; JURISDICTION AND VENUE: THIS WARRANTY, AND THE TERMS OF SALE OF THE PRODUCT, SHALL BE GOVERNED BY AND CONSTRUED IN ACCORDANCE WITH THE LAWS OF THE STATE OF FLORIDA WITHOUT REGARD TO PRINCIPLES OF CONFLICTS OF LAWS. THE PARTIES IRREVOCABLY SUBMIT AND CONSENT TO THE SOLE AND EXCLUSIVE JURISDICTION AND VENUE OF THE FLORIDA STATE COURTS IN AND FOR PALM BEACH COUNTY, FLORIDA, U.S.A. AND THE FEDERAL COURTS IN AND FOR THE SOUTHERN DISTRICT OF FLORIDA, U.S.A. THE PARTIES AGREE NOT TO RAISE THE DEFENSE OF FORUM NON CONVENIENS.

WARRANTY OUTSIDE THE U.S.A.

THE WARRANTY ON THIS PRODUCT SHALL COMPLY WITH APPLICABLE LAW WHEN SOLD BY AN AUTHORIZED DEALER TO A CONSUMER OUTSIDE OF THE UNITED STATES, INCLUDING EXTENDED WARRANTY COVERAGE PERIODS AS MAY BE LEGALLY REQUIRED IN SUCH JURISDICTIONS.

AUTHORIZED REPRESENTATIVE IN EUROPE (REGULATORY AFFAIRS ONLY):



Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

TEL: +31 (0) 70 345-8570 FAX: +31 70-346-7299

Actonis™

CUSTOMER SERVICE TOLL-FREE NUMBERS

United States and Canada 1-866-527-3726

United Kingdom 0808-134-9868

Australia 61 7 3198 2484

France 805-081-824

Spain 900-838-164

Mexico 01 800 953 0209

Italy 800 694 205

Germany 800 001 0163

**Authorized Representative in Europe
(Regulatory Affairs Only):**



Westervoortsedijk 60
6827 AT Arnhem
The Netherlands

Tel: +31 (0) 70 345-8570 Fax: +31 70 346-7299

**Australian Sponsor:
Emergo Australia**

Level 20, Tower II Darling Park,
201 Sussex Street Sydney,
NSW 2000 Australia



LEXINGTON INTERNATIONAL, LLC

1040 Holland Drive | Boca Raton, Florida | 33487 | USA
www.actonis.com