

### **How do I use the LED face mask?**

Step 1: Conduct a sensitivity test.

Step 2: Cleanse your skin and pat dry. Attach the mask, then adjust the strap around your head.

Step 3: Switch on, sit back and relax for your 10-minute treatment.

Step 4: Remove the mask from your head and wipe clean.

Step 5: Safely store in the protective storage bag.

### **Does it hurt to use?**

It should not feel uncomfortable.

### **How long will it take until I see results?**

You can expect to see visible results in just 4 weeks when used consistently 10 minutes a day, 3 to 5 times per week. Although effective on its own, you will achieve the best results when used with our CurrentBody Skin Hydrogel Face Masks or CurrentBody Skin Green Tea Serum.

**Can I use my normal moisturiser before I use the light therapy device?** Certain skincare and makeup can block the LED light from entering the skin. We recommend using the LED mask on clean, product free skin.

### **How many LEDs does the light therapy device have?**

There are 132 LED bulbs positioned in the face mask.

### **How long does the LED mask take to charge?**

6 to 8 hours for a full charge.

### **How should I clean the LED face mask after use?**

Wipe down the transparent side of the LED face mask with water-based wipes or a damp, lint-free cloth. Stubborn stains can be gently removed with 70% rubbing alcohol. We recommend cleaning your mask before every use to keep it fresh.

### **Does the CurrentBody Skin LED Mask have FDA approval?**

The CurrentBody LED Mask is FDA-cleared for at-home use.