

## **A710995 Hairmax LaserBand 82 ComfortFlex Hair Growth Device**

### **How to Use**

The Hairmax LaserBand 82 ComfortFlex should be used 3 times per week (every other day) regularly to achieve and maintain hair regrowth. Treatment time is as little as 90 seconds per treatment.

Fully charge the LaserBand for 4 to 6 hours before initial use (see User Manual for complete instructions). Use on damp or dry hair. Press the power button once to turn on. Using light pressure so the teeth can separate your hair, glide the LaserBand 82 into place at your hairline (this is position 1) and leave it there until it vibrates after 30 seconds. When you feel the vibration, glide the LaserBand back approximately 3 inches to position 2 (the middle portion of your scalp). Wait another 30 seconds until you feel the vibration again, then glide the LaserBand to the back of your scalp (position 3) and leave it for another 30 seconds. When your treatment is complete, turn the LaserBand off by pressing the power button once. The LaserBand will turn off automatically after 10 minutes.

For best results, charge after every use. Do not use the LaserBand 82 while it is plugged directly into the power cord.

To clean, make sure your device is unplugged and turned off. Wipe with a soft, slightly moistened cloth.

### **Where should I position the LaserBand to begin the treatment?**

The LaserBand should be placed at the front of your hairline to start your treatment. Cover your entire scalp during the recommended 90-second treatment time.

### **Do I have to keep using the Hairmax for the rest of my life?**

Hair regrowth with the LaserBand 82 is an ongoing process. For this reason, it is important that you use as directed on a consistent basis. Once you are satisfied with the results, you may reduce the frequency of treatments to once or twice per week for maintenance. It is not recommended that you discontinue the use of the Hairmax device altogether, as ongoing treatment is required to maintain benefits.

### **How often should the LaserBand be charged, and how should it be stored when not in use?**

For best results, charge after every use. The device can be stored on the charging cradle.

### **Can the LaserBand 82 be used by men and women?**

Yes, the LaserBand 82 can be used by men and women ages 18 and older.