

When & How



MINIMUM 2-3
TIMES PER WEEK



60 SECONDS PER
TARGET AREA



GLIDE TOOL IN ALL
DIRECTIONS



Face



Neck & Décolleté



Lips

1. Prior to GloPRO use, cleanse skin and prep with GloPRO Skin Prep Pads. The GloPRO Skin Prep Pads can be used on face and body to help condition the skin before microneedling.
2. Glide GloPRO on target areas, rolling in vertical, horizontal, and diagonal directions overlapping each area. Do NOT apply pressure; the roller will do the work for you! Let the weight of the tool guide your hand.
3. Dividing your face in four quadrants, roll over each quadrant for about 15 seconds. The GloPRO session for your entire face should last about 60 seconds. Treated areas may appear flushed after use. Clean MicroTip roller after every use by spraying with alcohol.
4. It is recommended for use two to three times a week; can be used daily.
5. Process is the same for face and/or body.

#GloTip: Tool can be used in the on or off mode. The LED light indicates that the product is turned on. Comes with two fully charged AAA batteries. Batteries are replaceable but not rechargeable.

Information provided by BeautyBio

*Do not use GloPRO on infected or irritated skin, open wounds or active breakouts.