

HOW TO USE

Emjoi® Micro-Pedi is designed for smoothing coarse skin on the feet and heels (fig. 6). Do not use this device anywhere else on the body.

Note: Only use Micro-Pedi roller on dry skin.

1. Wash and clean the skin and dry fully.

Note: Do not use on broken or bleeding skin.

2. Make sure the Micro-Pedi roller is placed into the unit. See section of **Changing the Micro-Pedi Roller** on how to remove and replace rollers.
3. Remove the protective cover from the unit and it is ready to use.



4. Turn on the unit by pressing the safety switch lock, in the center of the On/Off switch, and at the same time push the switch up (fig. 7). Gently rotate (back and forth or side to side) over the callous/corn/hard skin area for 2-3 seconds. **Do not press hard- gently glide the unit over desired area.** The unit will stop if it is pressed too hard on the skin. Stop and check if you have achieved the softness you desire. If not, go over the area for 2-3 more seconds and check again. Continue to use this method until the level of smoothness is obtained.

Warning: Stop using if skin becomes sore or inflamed.

Warning: Do not use on skin for more than 2-3 seconds at a time.

5. Turn the unit off (fig 8).
6. Rinse the skin or use a wet towel to get rid of the excess dead particles.
7. To clean the **Emjoi® Micro-Pedi**, please see the Cleaning section.
8. Keep the protective cover on your unit when it is not in use.

