


Using Your Elite 911

- 1** Choose a location for your charging cradle that's central to your home or living space and near a window. Plug in and charge the portable device for at least 3 hours. When it's fully charged, the light will turn solid red.



PLEASE NOTE: Do not plug the charging cradle power cord into an electrical outlet that is controlled by a light switch or into a power strip to avoid accidental power loss.

- 2** **Call 1 (800)797-1965**
to activate your device

or activate online at
[medicalguardian.com/retail-activation](https://www.medicalguardian.com/retail-activation)