

TROPHY SKIN

Diamond Tip Options



Standard Diamond Tip

Medium coarseness is suitable for face and body. Suction level on your device may be adjusted accordingly for more or less sensitive skin for an amazing microdermabrasion at home.



Fine Diamond Tip

The fine-grade tip is used to achieve spa-quality results at home without the discomfort of more abrasive tips.



Body Diamond Tip

Features a wider surface area with coarser diamond particles that was created to treat larger body parts such as hands, feet, legs, and chest.



Precision Diamond Tip

Perfect for around the eye area, this tip is used on the delicate, hard to reach, or angular areas of the face. Features smoother particles for a comfortable treatment in small spaces.